


# AKLARUS

DUAL WAVELENGTH  
PHOTOTHERAPY SYSTEM



The **AKLARUS Blue** and **Red** Combination System is an effective, drug-free alternative for treating acne and photodamaged skin. The non-invasive Aklarus treatment has been shown to significantly reduce inflammatory acne lesions. The narrow band blue light destroys the P. Acne bacteria while the Red light penetrates deeper stimulating cellular mechanisms responsible for tissue repair and regeneration.



Aklarus Anti-Aging **Infrared** and **Red** Combination System is used to treat wrinkles and photodamaged skin by stimulating collagen synthesis and increasing the amount of dermal elastic fibers.

## Aklarus Wavelength Specifications: **Blue, Red and Infrared**

Blue 420nm  
10W+/- 3W Total Output  
Standard dose 26J/ cm<sup>2</sup>  
in 20 Minutes

Red 628nm  
20W +/- 3W  
Standard dose 52J/cm<sup>2</sup>  
in 20 Minutes

Infrared 880nm  
16W +/- 3W  
Standard dose 42 J/cm<sup>2</sup>  
in 20 Minutes



# BEFORE & AFTER

AFTER 4 WEEKS / 8 TREATMENTS

These results have not been retouched in any way.

